

INSTRUCTIONS

FOR OPERATING VI-RAY-O HIGH FREQUENCY GENERATORS

READ CAREFULLY

First. Screw plug on end of cord into a convenient electric light socket. Then insert applicator into end of Generator (not too firmly.) Turn current on at lamp socket and Vi-Ray-O is ready for use.

Second. Turn adjusting screw at bottom to the right until you hear a humming sound and faintly see the Vi-Ray-O Rays in the applicator. Then place applicator upon that part of the body which is to be treated. To increase the current turn adjusting screw to the right. To decrease the current turn adjusting screw to the left.

Third. Start treatments with a weak current and gradually increase it.

NOTICE

Every Vi-Ray-O is thoroughly tested and adjusted before leaving the factory. Your Vi-Ray-O should work satisfactorily when you connect it to the lamp socket and turn the adjusting screw to the right.

In case it does not start, do not at once think the trouble is in the machine. It may be that the plug is not making good contact in your lamp socket. Be sure that the current is reaching the instrument.

This generator is fully guaranteed. It is not liable to get out of order, but in case it does, do not ship it to us until you have first written explaining as fully as possible how your Vi-Ray-O behaves. We may be able to tell you what your trouble is and thus save you needless express charges.

Under no circumstances take your Vi-Ray-O apart or permit anyone else to do so. It is at your own risk, if you do.

CAUTIONS

Do not use any hair tonic when treating the hair. Use only Vi-Ray-O High Frequency as a tonic.

Do not take a treatment in a bath tub.

Do not leave the instrument with current turned on—unscrew plug from lamp socket when through with treatment.

Do not bend cord sharply at point where it enters Vi-Ray-O. It will break in time, if you do.

Do not hang instrument up by the cord.

Do not turn adjusting screw with one hand while holding glass applicator with the other—a shock will be the result, if you do.

Eight to ten minutes should be the limit of any treatment. The instrument should be allowed to cool down before another application is made.

The vacuum in High Frequency applicators changes with prolonged use because some of the rarefied gas particles, disappearing under the action of Vi-Ray-O currents, become whiter and of a greenish yellow color. Any applicator which has reached this state can be re-exhausted by us at one-half of the original cost. Your Vi-Ray-O applicator is then as good as new.

All Vi-Ray-O applicators are made of the finest glass obtainable and should be handled with the greatest of care. Please note that it is impossible for us to guarantee applicators.

When ordering applicators state definitely number of applicator as given in our electrode circular which accompanies all machines. Add 10 cents to your remittance to cover parcel postage and insurance on your order.

WESTERN COIL & ELECTRICAL CO.

RACINE, WIS., U. S. A.

